

Voluntary Use of Respirators

This Notice is being provide to you in compliance with California Code of Regulations, Title 8, Appendix D to Section 5144: (Mandatory) Information for Employees Using Respirators When Not Required

Respirators are an effective method of protection against designated hazards when properly selected and worn. Respirator use is encouraged even when exposures are below the exposure limit, to provide an additional level of comfort and protection for workers. However, if a respirator is used improperly or not kept clean, the respirator itself can become a hazard to the worker. Sometimes, workers may wear respirators to avoid exposures to hazards, even if the amount of hazardous substance does not exceed the limits set by OSHA standards. If your employer provides respirators for your voluntary use, or if you provide your own respirator, you need to take certain precautions to be sure that the respirator itself does not present a hazard.

You should do the following:

1. Read and heed all instructions provided by the manufacturer on use, maintenance, cleaning and care, and warnings regarding the respirator's limitations.
2. Choose respirators certified for use to protect against the contaminant of concern. NIOSH, the National Institute for Occupational Safety and Health of the U.S. Department of Health and Human Services, certifies respirators. A label or statement of certification should appear on the respirator or respirator packaging. It will tell you what the respirator is designed for and how much it will protect you.
3. Do not wear your respirator into atmospheres containing contaminants for which your respirator is not designated to protect against. For example, a respirator designed to filter dust particles will not protect you against gases, vapors or very small solid particles of fumes or smoke.
4. Keep track of your respirator so that you do not mistakenly use someone else's respirator.

N95 Respirator

Summary

DMG

N95 Respirator

Foldable (Non-sterile)

BYD CARE 



- ⊕ Very close facial **fit**
- ⊕ $\geq 95\%$ filtrations efficiency against solid and liquid aerosols free of oil.



*Picture is for reference only, actual product may have differences.

Description

The N95 Respirator is designed to help provide respiratory protection for the user. This product has a filter efficiency of 95% or greater against solid and liquid aerosols free of oil. This product contains no components made from natural rubber latex.

Use instructions

- Suitable for protecting the mouth, nose, and chin to airborne particles, and block liquid aerosols, blood, body fluids, and secretions.
- For adult use only
- Do not wash.
- Store in a cool, dry, clean place away from fire and contamination.
- Store between -20°C and 38°C (-4°F and 100°F) with relative humidity below 80%.
- Expiration date: 2 years after production

3-Step check before use

1. Do not use if the packaging has been opened or damaged.
2. Check if the elastic bands are in good condition and are not damaged.
3. Check that the metallic strip is not broken.

Wearing instructions

1. Unfold the mask and hold it with both hands, with the metallic strip facing up.
2. Place elastic bands around the neck and head respectively. Put the mask against your face covering both nose and mouth.
3. Adjust the metallic strip over the bridge of the nose by using two fingers to press down until achieving a close fit.
4. Perform a fit test according to instructions.



Fit test

1. Place both hands over the mask.
2. Take a deep breath and hold your breath for a few seconds, making sure the mask collapses inward as you inhale.
3. Then exhale and hold your breath for a few seconds, making sure the mask bulges outward as you exhale.
4. If air leaks, reposition the mask, tighten the metallic strip and reposition the elastic bands for a better fit. Follow steps 1-3 again until a tight seal has been achieved.
5. Masks that have passed the fitting test are now safe to use.
6. Facial hair or other facial features may prevent the mask from achieving a proper fit. If a tight seal is not achievable, please DO NOT use this product.
7. If you experience difficulty breathing, dizziness, or other symptoms, leave the contaminated area immediately.